



# Oxygen Mask with Nebulizer

With Chamber and Tubing



Adult



Pediatric



Neonate



20180110



2018-01



2022-12

REGISTRO SANITARIO ECUADOR No. 348-DME-1014



YUYAO SHENGCHANG MEDICAL INSTRUMENT FACTORY

Yuyao, China

MADE IN CHINA



STERILE EO



0197

Direction for use:

1. Place prescribed medication in nebulizer jar by unscrewing lid.
2. Reattach lid.
3. Place the mask on the patient's face, centering both nose and mouth. The elastic strip should be below the ears and around the neck.
4. Gently pull the ends of the strap to adjust. Adjust metal strip at bridge of the nose to secure fit.
5. Attach nebulizer connection to face mask connection.

Caution:

- For single use only. Discard after use.
- Do not use if package is opened or damaged.
- Do not store at extreme temperature and humidity.  
Store in a cool and dry place.
- Ensure that all connections are secure and air circulates freely through the tubing.
- The tube, masks and nebulizers are latex free, only elastic belt is with latex.